





Welcome to the LabLink 360 interactive dashboard experience

Explore some of the exciting capabilities and conveniences in the all-new dashboard, part of the Premium Thermo Scientific™ LabLink™ 360 quality control software.

Click or tap each  below to learn more. Scroll down to explore the whole dashboard.


thermo scientific | LabLink 360 Quality Assurance Software
Q, Galaxy Mercury Lab (100002)
🔔

- Dashboard 
- Collapse 
- Review QC
- Levey-Jennings chart
- View peers
- Manage QC uploads
- Manage QC configurations
- Reports 
- Library 
- Settings
- Web Admin

Dashboard

Galaxy Mercury Lab (100002)

Welcome to LabLink 360 Quality Assurance Software
Clinical laboratory quality control software to simplify your QC analysis


QC review 

5


Review 1 →


Approve 4 →


Awaiting monthly approval 0 →

QC alerts 

3

Warnings  1 →

Exclusions  2 →

Missed QC 


0

Instruments

See more details

Select a date:

Last 24 hours 7 days 14 days 30 days




Auto-communicator connection to LabLink 360

System connected

Last time connected: 09/26/2025 at 04:49:10 am

AC software → LabLink 360

Quality Control Performance 

Laboratory, Instrument, and Assay data

Laboratory

1179
QC results

30 days

Total QC results

- Within control limits 1025 87%
- Warnings 53 4%
- Exclusions 101 9%

30 days

QC results

Comments 89

Events 4

General

Comments 3


Events 10

Last six months

Warnings and exclusions

Month	Count	% Change
September	583	
August	550	-31.34%
July	550	14.43%
June	700	11.29%
May	629	-3.53%
April	652	-15.32%

Last Updated 12:04:13 pm. [Update](#)

Instruments 


Total: 4 instruments

ROCHE COBAS 4000, Cobas 1

Siemens Atellica CH, Atellica 1

Siemens Atellica CH, Atellica 3

Number of QC results

Assays 

Total: 8 assays (6 products)

Bilirubin, Di...

Bilirubin, Tot...

Ammonia


Albumin

Acetaminophen

Ethanol (EtOH)


Ammonia

C-Reactive Pro...

QC uploads 

View the import status for recently uploaded files


File name	Status	Date	Time	Upload routine
TA289435.CSV	Parser issue	09/12/2025	11:36:07 am	Auto-Communicator
A004562.CSV	Imported	09/12/2025	11:35:37 am	Auto-Communicator
DL0289200.CSV	Imported	09/12/2025	11:34:28 am	Auto-Communicator
LS292800.CSV	Imported	09/12/2025	11:21:01 am	Auto-Communicator
LS294345.CSV	Imported	09/12/2025	11:20:20 am	Auto-Communicator
TA267492.CSV	Mapping issue	09/12/2025	11:13:17 am	Auto-Communicator

Documents 

6 new documents

Updates to subscribed products and instruments


- MAS™ Alcohol and Ammonia Controls
- MAS™ Omni-CORE Controls
- MAS™ Omni-IMMUNE™ Immunoassay
- MAS™ Omni-CORE Controls
- MAS™ Alcohol and Ammonia Controls
- MAS™ Alcohol and Ammonia Controls
- MAS™ Alcohol and Ammonia Controls
- MAS™ Alcohol and Ammonia Controls
- MAS™ Alcohol and Ammonia Controls

Product lots 

18 lots about to expire

Expiring within 90 days

Lot number	Product	Expiration date
BC25041A	MAS Bilirubin	11/23/2025
BC25042A	MAS Bilirubin	11/23/2025
BC25043A	MAS Bilirubin	11/23/2025
CHU25121A	MAS ChemTRAK-H Unassayed	11/23/2025
CHU25122A	MAS ChemTRAK-H Unassayed	11/23/2025
CHU25123A	MAS ChemTRAK-H Unassayed	11/23/2025
UA28021A	MAS Liquimmune	11/23/2025
IM25081A	MAS Immunology	11/23/2025


Reports 

Newest reports

Go to reports activity

Report Type	Created date and time
Levey-Jennings	09/26/2025 at 06:09:58 am
Method Performance Sigma	09/25/2025 at 07:06:32 pm
Affiliated Lab	09/17/2025 at 07:30:46 am
Statistical Summary	09/15/2025 at 11:43:54 pm
Monthly Summary	09/13/2025 at 12:41:55 am
Exception Notes	09/13/2025 at 12:31:02 am
Individual Lab	09/13/2025 at 12:23:20 am

Peer Comparison: A Pragmatic Approach by Sten Westgard



Labs need Peer Group monitoring due to infrequent PT/EOA. But reports overwhelm with stats & graphs. Key questions: What numbers matter? Is SDI good or just okay? How many graphs? Does group size matter for effective peer comparison?

[Register to watch it on-demand](#)

Recommended products

MAS™ Diabetes™ Max Quality Controls

MAS™ Diabetes™ Max

MAS™ Omni-CORE™ Controls

MAS™ Omni-CORE™

MAS™ Omni-CARDIO™ Control

MAS™ Omni-CARDIO™ Control

MAS™ ChemTRAK™ H General Chemistry Controls

MAS™ ChemTRAK™ H

© 2025 Thermo Fisher Scientific. All rights reserved. | [Privacy policy](#) | [Terms of use](#) | [Terms of service](#) Version 1.0.209 / validation

Accelerate your QC review by upgrading to LabLink 360 today!

For more information, please contact your Thermo Scientific sales representative.

Learn more at thermofisher.com/lablink360

thermo scientific

Not all products are CE marked or have FDA clearance for sale in the U.S. Availability of products in each country depends on local regulatory marketing authorization status. © 2025 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. CDD-FR-MTL-0687